

# graphical\_user\_interface

*Pamgene Support*

*6/19/2019*

## Introduction

A user requirement for a pamservice tool to be used with

- existing ps12 instruments
- new ps12 instruments


The requirements cover the graphical user interface (GUI) for servicing the ps12.

## Graphical User Interface

The requirements for a the GUI is outlined below and each tab is dealt with in each section.

### main tab

Used to decide which board to connect to. This tab is removed in the new pamservice tool.



### parameters tab

used to get and set the values of the parameters

Main	Parameters	Disposable	Pumping	Position	Camera	Dispense	Input
------	------------	------------	---------	----------	--------	----------	-------

Controllable Settings	Sensor Calibration	Logical Positions
-----------------------	--------------------	-------------------

eEncFilterB

rPosFilter1	rPosFilter2	rPosFilter3
<input type="text" value="0,00"/>	<input type="text" value="0,00"/>	<input type="text" value="0,00"/>

eEncFocusB

rPosForFilter1	rPosForFilter2	rPosForFilter3	rOffsetForDisp1	rOffsetForDisp2	rOffsetForDisp3
<input type="text" value="0,00"/>	<input type="text" value="0,00"/>	<input type="text" value="0,00"/>	<input type="text" value="0,00"/>	<input type="text" value="0,00"/>	<input type="text" value="0,00"/>

eEncTransXB

xyLoadPosition	<input type="text" value="0,00"/>
xyIncubationPosition	<input type="text" value="0,00"/>
xyFrontPanelPosition	<input type="text" value="0,00"/>
xyReadPosition	<input type="text" value="0,00"/>
xyAspiratePosition	<input type="text" value="0,00"/>
xyDispensePosition	<input type="text" value="0,00"/>
xyOffset_Well	<input type="text" value="0,00"/>
xyOffset_Disposable	<input type="text" value="0,00"/>
xyOffset_FocusGlass	<input type="text" value="0,00"/>
xyOffset_DispenseHead1	<input type="text" value="0,00"/>
xyOffset_DispenseHead2	<input type="text" value="0,00"/>
xyOffset_DispenseHead3	<input type="text" value="0,00"/>
xyOffset_DispenseHead4	<input type="text" value="0,00"/>
xyOffset_Prime1	<input type="text" value="0,00"/>
xyOffset_Prime2	<input type="text" value="0,00"/>

eEncTransYB

xyLoadPosition	<input type="text" value="0,00"/>
xyIncubationPosition	<input type="text" value="0,00"/>
xyFrontPanelPosition	<input type="text" value="0,00"/>
xyReadPosition	<input type="text" value="0,00"/>
xyAspiratePosition	<input type="text" value="0,00"/>
xyDispensePosition	<input type="text" value="0,00"/>
xyOffset_Well	<input type="text" value="0,00"/>
xyOffset_Disposable	<input type="text" value="0,00"/>
xyOffset_FocusGlass	<input type="text" value="0,00"/>
xyOffset_DispenseHead1	<input type="text" value="0,00"/>
xyOffset_DispenseHead2	<input type="text" value="0,00"/>
xyOffset_DispenseHead3	<input type="text" value="0,00"/>
xyOffset_DispenseHead4	<input type="text" value="0,00"/>
xyOffset_Prime1	<input type="text" value="0,00"/>
xyOffset_Prime2	<input type="text" value="0,00"/>

eLEDUnit

Intensity_LED1	<input type="text" value="0"/>
Intensity_LED2	<input type="text" value="0"/>
Intensity_LED3	<input type="text" value="0"/>


Read from file	Write to file	Flash Rest	Flash Sens	Read All
----------------	---------------	------------	------------	----------

Controllable Settings		Sensor Calibration			Logical Positions		
Name		eSensorObj	rSysOffset	rSysGain	rSysGain2	rSensOffset	rSensGain
0	eTempSensDispense	108	9,997600E+1	3,050500E-3	1,379500E-13	-2,448300E+2	2,3
	eTempSensBottomLeft	50	9,999700E+1	3,051600E-3	1,825800E-11	-2,456000E+2	2,3
	eTempSensBottomMiddle	52	9,999400E+1	3,050100E-3	8,503900E-11	-2,469600E+2	2,3
	eTempSensBottomRight	54	1,000200E+2	3,050500E-3	1,245400E-10	-2,475900E+2	2,3
	eTempSensCoverLeft	49	1,000070E+2	3,050800E-3	5,956300E-11	-2,435200E+2	2,3
	eTempSensCoverMiddle	51	9,999100E+1	3,051400E-3	1,447000E-12	-2,482500E+2	2,4
	eTempSensCoverRight	53	1,000080E+2	3,050900E-3	8,379600E-11	-2,473300E+2	2,3
	eTempSensCalibration	55	1,000000E+2	3,100000E-3	0,000000E+0	-2,454500E+2	2,3
	ePressSensUnderPress	2	-4,236000E-1	1,100000E+0	0,000000E+0	0,000000E+0	-2,2
	ePressSensAspirate	109	-2,091000E-1	1,000000E+0	0,000000E+0	0,000000E+0	-2,2
	ePressSensOverPress	1	-1,182000E-1	1,050000E+0	0,000000E+0	0,000000E+0	2,2
	ePressSensWell1	56	-2,245500E+0	1,000000E+0	0,000000E+0	0,000000E+0	4,4
	ePressSensWell2	57	-2,268200E+0	1,000000E+0	0,000000E+0	0,000000E+0	4,4
	ePressSensWell3	58	-2,222700E+0	1,000000E+0	0,000000E+0	0,000000E+0	4,4

Controllable Settings		Sensor Calibration			Logical Positions		
Name		eCtrlObj	eSensObj	eCurrSensObj	uWaitReadyTimeout	iDefaultPwm	rMaxP
0	eCtrlDispenseHeadHeater	138	108	110	900000	0	90,00
	eCtrlHtrCoverLeft	101	49	71	900000	0	88,44
	eCtrlHtrCoverMiddle	103	51	72	900000	0	88,45
	eCtrlHtrCoverRight	105	53	73	900000	0	88,48
	eCtrlHtrBottomLeft	102	50	68	900000	0	88,45
	eCtrlHtrBottomMiddle	104	52	69	900000	0	88,38
	eCtrlHtrBottomRight	106	54	70	900000	0	88,44
	eCtrlUnderPressureSupplier	43	2	0	10000	0	0,02
	eCtrlOverPressureSupplier	42	1	0	10000	0	0,60
	eCtrlAspirateHeadPositioner	137	128	111	80000	30	2,00
	eCtrlDispenseHeadPositioner1	139	112	0	30000	100	100,00
	eCtrlDispenseHeadPositioner2	140	114	0	30000	100	100,00
	eCtrlDispenseHeadPositioner3	141	116	0	30000	100	100,00
	eCtrlDispenseHeadPositioner4	142	118	0	30000	100	100,00
	eCtrlFilterPositioner	47	16	6	30000	0	10,00


## disposable tab

used to set and get the disposable parameters



## position tab

used to set and get of the position parameters




## pumping tab

used to set and get of the pumping parameters

Main	Parameters	Disposable	Pumping
<div><div>eOverPressSup</div><div>Target0.20</div><div>Actual0.20</div><div>HeatRamp0.00</div><div>E R Err H</div><div><div></div><div></div><div></div><div></div></div></div> <div><div>eUnderPressSup</div><div>Target0.09</div><div>Actual0.01</div><div>HeatRamp0.00</div><div>E R Err H</div><div><div></div><div></div><div></div><div></div></div></div>			


## camera tab

used to control the camera



## dispense tab

Used to have setting and getting of dispense.





## input/outputs tab

Used to get and set the values across the three control boards.

- iccb (incubate board)
- sycb (system board)
- wsb (wash board)

## iccb (incubate board)

The following four tabs are moved into one incube board tab.



WSCB Analog In		ICCB Analog Out	
Digital In		ICCB Digital Out	
eFanCover	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eFanBottom	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell1	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell2	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell3	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell4	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell5	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell6	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell7	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell8	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell9	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell10	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell11	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveWell12	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValvePrsSelection	<input checked="" type="checkbox"/>		<input type="checkbox"/>
eValveVacSelection	<input checked="" type="checkbox"/>		<input type="checkbox"/>

Select all values button



log In		WSCB Analog In		ICCB Analog Out		SYCB Analog Out		WSCB Analog Out	
al In		SYCB Digital In		WSCB Digital In		ICCB Digital Out		SYCB Digital Out	
		WSCB Digital Out		ICCB Analog In					
eTempSensCoverLeft	0,00			eCurrSensBottomMiddle	0,00				
eTempSensBottomLeft	0,00			eCurrSensBottomRight	0,00				
eTempSensCoverMiddle	0,00			eCurrSensCoverLeft	0,00				
eTempSensBottomMiddle	0,00			eCurrSensCoverMiddle	0,00				
eTempSensCoverRight	0,00			eCurrSensCoverRight	0,00				
eTempSensBottomRight	0,00								
eTempSensCalibration	0,00								
ePressSensWell1	0,00								
ePressSensWell2	0,00								
ePressSensWell3	0,00								
ePressSensWell4	0,00								
ePressSensWell5	0,00								
ePressSensWell6	0,00								
ePressSensWell7	0,00								
ePressSensWell8	0,00								
ePressSensWell9	0,00								
ePressSensWell10	0,00								
ePressSensWell11	0,00								
ePressSensWell12	0,00								
eCurrSensBottomLeft	0,00								

sycb (system board)

The following four tabs in one incube board tab.

ICCB Digital In		SYCB Digital In		WSCB Digital In	
SYCB Analog In		WSCB Analog In			
ePressSensOverPress	0,00				
ePressSensUnderPress	0,00				
eCurrSensTransX	0,00				
eCurrSensTransY	0,00				
eCurrSensFocus	0,00				
eCurrSensFilter	0,00				

ICCB Digital In	SYCB Digital In	WSCB Digital In	ICCB Digital Out	SYCB Digital Out	WSCB Digital Out
SYCB Analog In	WSCB Analog In	ICCB Analog Out	SYCB Analog Out		


  

ePWMSpare1	0,00		0,00
eMotTransporterX	0,00		0,00
eMotTransporterY	0,00		0,00
eMotFocus	0,00		0,00
eMotFilter	0,00		0,00

SYCB Analog In	WSCB Analog In	ICCB Analog In
ICCB Digital In	SYCB Digital In	WSCB Digital In


  

eEncTransXA	
eEncTransXB	
eEncTransXI	
eEncTransYA	
eEncTransYB	
eEncTransYI	
eEncFocusA	
eEncFocusB	
eEncFocusI	
eEncFilterA	
eEncFilterB	
eEncFilterI	
eSensTransXHome	
eSensTransYHome	
eSensFocusHome	
eSensFilterHome	
eSensOISpare1	
eSensSysPowerOK	
eDiagSYCBPOError	



## wsb (wash board)

The following four tabs are moved into one wash board tab.



SYCB Analog In	WSCB Analog In	ICCB
eTempSensDispense	0,00	
ePressSensAspirate	0,00	
eCurrSensDispenseHeater	0,00	
eCurrSensAspirateMotor	0,00	

ICCB Digital In	SYCB Digital In	WSCB Digital In	ICCB Digital Out	SYCB Digital Out
SYCB Analog In	WSCB Analog In	ICCB Analog Out	SYCB Analog Out	
ePWMSpare1	0,00	$\frac{f}{T}$	0,00	
eMotTransporterX	0,00	$\frac{f}{T}$	0,00	
eMotTransporterY	0,00	$\frac{f}{T}$	0,00	
eMotFocus	0,00	$\frac{f}{T}$	0,00	
eMotFilter	0,00	$\frac{f}{T}$	0,00	

WSCB Digital In	ICCB Digital Out	SYCB Digital Out	WSCB Digital Out
ePmpAspirate	<input checked="" type="checkbox"/>		

error tab

